



BODYSTEEL & SILVER CO.,LTD
63 Moo 5, Soi Wat Bang-Co Yaek Rong-Foam, Kanchanaphisek Road
Bang Ma-Nang, Bang Yai, NONTABURI 11140 THAILAND

☎ 66 2-1505999 📠 66-2-1505990

sales@steelandsilver.com

EO STERILIZATION Quality Control Report

Date : **15-FEB-2024**

LOT : 0926

STERILIZATION DATE : 07-FEB-2024

EXPIRE DATE : 07-FEB-2029

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.

Result
PASS



Sample Fail



Yellow

Sample New
Indicator



Pink

Quality Assurance By:

Date: **15-FEB-2024**



STERIL 45 C DAT : 0926		07/02/2024 14:43		STEP TIME		TEMP(C)		PRES(KPa)		HUMI(20H)	
PH	14:43	28.4	00	00	89.8						
PV	16:15	45.0	00	00	83.3						
HU	16:28	44.6	-85	-85	75.9						
OS	16:38	44.9	-84	-84	78.5						
ST	16:40	45.4	-50	-50	NONE						
ST	17:10	44.9	-48	-48	NONE						
ST	17:40	44.6	-48	-48	NONE						
ST	18:10	45.3	-47	-47	NONE						
ST	18:40	45.0	-47	-47	NONE						
ST	19:10	44.6	-47	-47	NONE						
ST	19:40	45.3	-47	-47	NONE						
ST	20:10	45.4	-47	-47	NONE						
ST	20:40	44.8	-47	-47	NONE						
ST	21:10	44.7	-47	-47	NONE						
AE	21:40	45.4	-46	-46	NONE						
AE	22:39	45.1	-88	-88	NONE						
AE	23:40	44.8	-87	-87	NONE						
AE	00:41	45.5	-88	-88	NONE						
AE	01:41	44.8	-87	-87	NONE						
AE	02:42	45.6	-88	-88	NONE						
AE	03:42	44.8	-88	-88	NONE						
AE	04:43	45.5	-88	-88	NONE						
AE	05:43	44.8	-88	-88	NONE						
AE	06:43	45.5	-88	-88	NONE						
AE	09:27	35.0	-88	-88	NONE						
AE	09:27	35.8	-88	-88	NONE						
FI	09:37	37.7	-84	-84	NONE						
OK BY :											
APP BY :											
RESULT :											

Protocol: