



BODYSTEEL & SILVER CO.,LTD
63 Moo 5, Soi Wat Bang-Co Yaek Rong-Foam, Kanchanaphisek Road
Bang Ma-Nang, Bang Yai, NONTABURI 11140 THAILAND

☎ 66 2-1505999 📠 66-2-1505990

sales@steelsilver.com

EO STERILIZATION Quality Control Report

Date : **15-FEB-2024**

LOT : 0927

STERILIZATION DATE : 08-FEB-2024

EXPIRE DATE : 08-FEB-2029

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.

Result
PASS



Sample Fail



Yellow

Sample New Indicator



Pink

Quality Assurance By:

Date: **15-FEB-2024**



-STERI #45 G DAY 0927		08/02/2024 16:03		STEP TIME		TEMP (C)		PRESS (KPa)		HUMID (GRH)	
PH		16:03		23.9		00		87.6			
17:26	45.0	00	65.4	17:39	44.8	-85	77.5	17:49	44.7	-84	80.5
17:52	45.2	-50	NONE	18:22	45.0	-49	NONE	18:52	44.6	-49	NONE
19:22	45.1	-49	NONE	19:52	45.5	-48	NONE	20:22	44.7	-48	NONE
20:52	44.7	-48	NONE	21:22	45.5	-48	NONE	21:52	45.3	-48	NONE
22:52	44.6	-48	NONE	23:51	44.9	-87	NONE	00:53	45.7	-88	NONE
01:53	44.8	-88	NONE	02:54	45.7	-88	NONE	03:55	44.8	-88	NONE
04:55	45.6	-88	NONE	05:55	44.8	-87	NONE	06:55	45.6	-88	NONE
07:54	45.0	-88	NONE	09:53	40.0	-82	NONE	09:55	39.9	-87	NONE
CHK BY :		APP BY :		RESULT :							

Protocol: