



BODYSTEEL & SILVER CO.,LTD
63 Moo 5, Soi Wat Bang-Co Yaek Rong-Foam, Kanchanaphisek Road
Bang Ma-Nang, Bang Yai, NONTHABURI 11140 THAILAND

☎ 66 2-1505999 📠 66-2-1505990

sales@steelsilver.com

EO STERILIZATION Quality Control Report

Date : 22-FEB-2024

LOT : 0929

STERILIZATION DATE : 19-FEB-2024

EXPIRE DATE : 19-FEB-2029

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.

Result
PASS



Sample Fail



Yellow

Sample New
Indicator



Pink

Quality Assurance By:

Date: 22-FEB-2024



STERI : 45 C. DAT 19/02/2024 16:02	STEP TIME	TEMP (C)	PRES (KPa)	HM1 (20H)	HM2 (20H)
PH 16:02	28.7	00	88.9		
FW 17:05	45.0	00	86.3		
HU 17:48	45.2	-85	76.2		
GS 17:58	44.6	-84	79.8		
ST 18:00	45.3	-50	NONE		
ST 18:30	45.3	-48	NONE		
ST 19:00	44.6	-48	NONE		
ST 19:30	45.1	-47	NONE		
ST 20:00	45.3	-47	NONE		
ST 20:30	44.7	-47	NONE		
ST 21:00	44.9	-47	NONE		
ST 21:30	45.5	-47	NONE		
ST 22:00	45.1	-47	NONE		
ST 22:30	44.7	-47	NONE		
AE 23:00	45.3	-47	NONE		
AE 23:30	45.6	-88	NONE		
AE 01:00	45.8	-88	NONE		
AE 02:00	44.8	-88	NONE		
AE 03:01	45.7	-87	NONE		
AE 04:01	44.8	-88	NONE		
AE 05:02	45.6	-88	NONE		
AE 06:02	44.8	-88	NONE		
AE 07:02	45.7	-88	NONE		
AE 08:01	44.8	-88	NONE		
AE 09:03	00.0	-83	NONE		
AE 09:05	43.8	-87	NONE		
FI 09:14	44.6	-80	NONE		
OK BY :					
APP BY :					
RESULT :					

Protocol: