



BODYSTEEL & SILVER CO.,LTD
63 Moo 5, Soi Wat Bang-Co Yaek Rong-Foam, Kanchanaphisek Road
Bang Ma-Nang, Bang Yai, NONTHABURI 11140 THAILAND

☎ 66 2-1505999 📠 66-2-1505990

sales@steelsilver.com

EO STERILIZATION Quality Control Report

Date : 29-MAR-2024

LOT : 0937

STERILIZATION DATE : 26-MAR-2024

EXPIRE DATE : 26-MAR-2029

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.

Result
PASS



Sample Fail



Yellow


Sample New Indicator



Pink

Quality Assurance By:

Date: 29-MAR-2024



LOT 937									
Parameters									
Start	End	Temp	Time	Pressure	Humidity	Flow	Leak	Alarm	Notes
31.0	31.1	-0.2	98.9	14.18					
Water	31.1	-0.2	74.7	14.18					
Vacuum	31.1	-5.3	74.1	14.19					
Preheat	35.0	-3.8	74.3	14.26					
Heat	38.7	-2.6	74.1	14.36					
Vacuum	37.3	-60.0	74.0	14.41					
Keep	37.9	-59.1	74.5	14.49					
Intake	37.8	-5.0	74.9	14.52					
Vacuum	37.1	-60.1	74.7	14.57					
Humid	37.1	-60.1	74.7	14.57					
Humid	36.9	-59.7	74.7	15.27					
Dry	36.7	-57.7	74.0	15.36					
Post	36.9	-51.5	74.7	15.36					
Diff	38.3	-24.6	74.9	15.46					
Peak	50.1	-16.0	74.8	15.20					
Start	50.0	-16.0	74.6	16.20					
Start	51.4	-13.7	74.6	16.30					
Start	50.5	-11.7	74.2	16.40					
Start	50.4	-10.1	74.0	16.50					
Start	50.3	-9.7	74.6	17.00					
Start	50.2	-9.3	74.8	17.10					
Start	50.2	-9.0	74.4	17.20					
Start	50.3	-8.6	74.4	17.30					
Start	50.1	-8.5	74.4	17.40					
Start	50.4	-8.0	74.0	17.50					
Start	50.3	-7.8	74.8	18.00					
Start	50.1	-7.5	74.4	18.10					
Start	50.3	-7.3	74.4	18.20					
Start	50.3	-7.0	74.2	18.30					
Start	50.1	-6.8	74.0	18.40					
Start	50.3	-6.6	74.2	18.50					
Start	50.1	-6.4	74.2	19.00					
Start	50.1	-6.2	74.4	19.10					
Start	50.1	-6.0	74.1	19.20					
Start	50.2	-5.7	74.2	19.30					
Start	50.3	-5.4	74.8	19.40					
Start	50.1	-5.1	74.0	19.50					
Start	50.2	-4.8	74.0	20.00					
Start	50.1	-4.5	74.4	20.10					
Start	50.3	-4.1	74.0	20.20					
Start	50.4	-3.5	74.0	20.30					
Start	50.3	-3.1	74.2	20.40					
Start	50.3	-2.7	74.6	20.50					
Start	50.0	-2.4	74.6	21.00					
Start	50.2	-2.0	74.8	21.10					
Start	50.0	-1.7	74.8	21.20					
Start	50.3	-1.5	74.0	21.30					
Start	50.3	-1.3	74.0	21.40					
Start	50.2	-1.1	74.4	21.50					
Start	50.4	-1.0	74.2	22.00					
Start	50.4	-0.9	74.4	22.10					
Start	50.2	-0.9	74.0	22.20					
Clean	51.3	-17.1	74.9	22.53					
Intake	50.3	-16.0	74.2	23.24					
Diff	51.1	-16.4	74.2	23.56					
4th	51.4	-16.9	74.6	00.27					
5th	51.9	-16.9	74.7	00.50					
Intake	51.3	-10.1	74.1	00.59					
End	51.3	-10.1	74.1	00.59					

Protocol: