



Date : 03-Oct-2025

LOT : 1155

STERILIZATION DATE : 20-Sep-2025

EXPIRE DATE : 20-Sep-2029

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.



Sample Fail



Yellow

Sample New
Indicator



Pink

Quality Assurance By:

Date: 03-Oct-2025



Lot 1065									
Interno use only: send medical exp.									
Ethyline oxide sterilize									
Start:	50 °C								
Stop:	55 °C								
Start time:	6h	0min							
Client time:	6 p/c								
Client p:	0 kPa								
Client time:	15 min								
T °C P kPa H kPa time									
Start:	24.8	-0.6	95.9	15.40					
Water:	24.8	-0.5	74.7	15.41					
Vacu:	24.8	-5.1	74.8	15.41					
Preh:	35.0	-1.7	74.2	15.56					
KHeat:	35.1	-1.6	74.4	16.06					
Vacu:	33.8	-8.0	74.3	16.12					
Washp:	35.9	-57.9	74.3	16.20					
Start:	35.0	-4.9	74.8	16.23					
Vacu:	38.3	-80.1	74.7	16.28					
Washp:	38.3	-80.1	74.8	16.28					
Start:	37.3	-67.8	74.7	16.36					
Washp:	37.3	-67.8	74.7	16.36					
Start:	37.3	-48.4	74.2	17.06					
Diffu:	37.4	-24.4	74.2	17.16					
Start:	45.9	-12.4	74.6	17.53					
Washp:	45.9	-12.4	74.6	17.53					
Start:	51.4	-9.1	74.4	18.03					
Start:	50.9	-7.8	74.4	18.03					
Start:	50.3	-6.8	74.4	18.23					
Start:	50.3	-5.2	74.8	18.43					
Start:	50.2	-5.2	74.8	18.43					
Start:	50.3	-4.3	74.8	18.53					
Start:	50.3	-4.3	74.8	18.53					
Start:	50.3	-2.9	74.6	19.13					
Start:	50.1	-2.0	74.6	19.23					
Start:	50.3	-1.4	74.4	19.33					
Start:	50.1	-1.0	74.4	19.43					
Start:	50.3	-0.8	74.4	19.53					
Start:	50.3	-0.5	74.8	20.03					
Start:	50.1	-0.5	74.8	20.13					
Start:	50.4	-0.4	74.0	20.23					
Start:	50.0	-0.4	74.6	20.33					
Start:	50.3	-0.4	74.6	20.43					
Start:	50.1	-0.3	74.6	20.53					
Start:	50.1	-0.3	74.6	20.53					
Start:	50.3	-0.4	74.4	21.13					
Start:	50.1	-0.4	74.4	21.23					
Start:	50.1	-0.4	74.2	21.33					
Start:	50.3	-0.4	74.0	21.43					
Start:	50.2	-0.5	74.6	21.53					
Start:	50.2	-0.5	74.2	22.03					
Start:	50.1	-0.5	74.4	22.13					
Start:	50.1	-0.5	74.4	22.13					
Start:	50.4	-0.6	74.6	22.23					
Start:	50.3	-0.6	74.4	22.33					
Start:	50.3	-0.6	74.4	22.43					
Start:	50.3	-0.6	74.4	22.43					
Start:	50.1	-0.6	74.2	23.03					
Start:	50.1	-0.6	74.2	23.03					
Start:	50.3	-0.6	74.4	23.13					
Start:	50.3	-0.6	74.6	23.23					
Start:	50.0	-0.6	74.0	23.33					
Start:	50.3	-0.6	74.2	23.43					
Start:	50.2	-0.6	74.3	23.53					
Clean Process									
1th:	50.7	-15.9	74.6	00:19					
2th:	50.9	-15.8	74.6	01:46					
3th:	50.9	-15.6	74.6	01:12					
4th:	50.9	-15.6	74.4	01:58					
5th:	50.7	-15.6	74.4	01:58					
6th:	50.4	-15.8	74.4	02:32					
7th:	50.5	-8.9	74.3	02:33					
End:	50.5	-8.9	74.3	02:33					
Sterilization Info Completed									
Steril Time:	106520m								
Issued:	00217								
Date:	2025-09-20								
Number:	50000450051								
Operator:									

Protocol: