



BODYSTEEL & SILVER CO.,LTD
63 Moo 5, Soi Wat Bang-Co Yaek Rong-Foam, Kanchanaphisek Road
Bang Ma-Nang, Bang Yai, NONTABURI 11140 THAILAND

☎ 66 2-1505999 📠 66-2-1505990

sales@steelsilver.com

EO STERILIZATION Quality Control Report

Date : 7-Aug-2026

LOT : 1261

STERILIZATION DATE : 28-July-2026

EXPIRE DATE : 28-July-2031

QC REPORT :

PASS as the test result has been shown as pictures below.

BIOLOGICAL INDICATORS :

The EZTest indicator is placed in a suitable test pack which is representative of the load. The load is then processed as usual.

After sterilization the chemical indicator ring on the label changes from blue to brownish when exposed to Ethylene Oxide (EO gas). This distinguishes exposed from unexposed units (it does not indicate acceptable sterilization).

After the incubation period of 48 hours, the appearance of a yellow color indicates bacterial growth. No color change indicates adequate sterilization.

Result
PASS



Sample Fail



Yellow

Sample New Indicator



Pink

Quality Assurance By:

Date: 7-Aug-2026

W. Tamee



Welcome use sampling medical equ Sterilize on the sterilite		LOT 1261		T C P kPa H % time	
Start:	25.2	-0.6	3.2	16:44	
Water:	25.2	-0.5	35.0	16:45	
Vacuum:	25.1	-5.1	35.6	16:45	
Preheat:	35.0	-3.3	35.0	17:00	
Heat:	35.2	-2.6	35.0	17:10	
Vacuum:	35.6	-6.0	35.0	17:16	
Keep:	35.6	-5.8	35.2	17:24	
Inta:	40.0	-4.9	35.3	17:31	
Vacuum:	37.0	-5.0	35.6	17:46	
Humid:	42.1	-5.6	35.2	17:46	
Post:	39.3	-5.7	35.4	18:15	
Diffu:	38.9	-4.9	35.1	18:23	
Heat:	50.1	-11.8	35.3	18:38	
Steri:	50.0	-11.8	35.5	19:17	
Steri:	51.3	-9.6	35.9	19:27	
Steri:	50.6	-8.7	35.3	19:37	
Steri:	50.8	-8.0	35.3	19:47	
Steri:	50.4	-7.5	35.3	19:57	
Steri:	50.3	-7.0	35.1	20:07	
Steri:	50.1	-6.3	35.7	20:17	
Steri:	50.3	-6.1	35.5	20:27	
Steri:	50.0	-5.9	35.9	20:47	
Steri:	50.4	-5.6	35.1	20:57	
Steri:	50.0	-5.4	35.2	21:07	
Steri:	50.4	-5.2	35.7	21:17	
Steri:	50.2	-5.1	35.5	21:27	
Steri:	50.1	-4.8	35.8	21:37	
Steri:	50.2	-4.7	35.9	21:47	
Steri:	50.4	-4.6	35.5	21:57	
Steri:	50.0	-4.5	35.5	22:07	
Steri:	50.1	-4.4	35.3	22:17	
Steri:	50.4	-4.3	35.5	22:27	
Steri:	50.2	-4.2	35.1	22:47	
Steri:	50.3	-4.1	35.1	23:07	
Steri:	50.2	-4.0	35.9	23:17	
Steri:	50.4	-3.9	35.9	23:27	
Steri:	50.0	-3.9	35.7	23:37	
Steri:	50.4	-3.9	35.9	23:47	
Steri:	50.3	-3.8	35.3	00:07	
Steri:	50.3	-3.7	35.9	00:27	
Steri:	50.2	-3.7	35.7	00:37	
Steri:	50.4	-3.7	35.9	00:47	
Steri:	50.1	-3.7	35.1	01:07	
Steri:	50.4	-3.6	35.0	01:17	
Clean:	50.0	-3.6	35.0	01:17	
11h22am	50.0	-3.6	35.7	01:43	
2th	50.4	-40.7	35.1	02:13	
3th	50.8	-40.3	35.3	02:44	
4th	51.2	-40.1	35.2	03:15	
5th	51.4	-40.5	35.6	03:45	
6th	51.5	-40.3	35.0	04:16	
Inta:	51.4	-35.9	35.8	04:17	
End:	51.4	-35.9	35.8	04:17	
Sterilization Completed					
Steril Time: 11h22am					
Usecount: 00249					
Date: 2026-7-28					
Number: 5000045005					
Operator:					

Protocol: